

Instructions for Use

1. Remove the battery compartment cover by lightly depressing the cover with your thumb and slide it in the direction of the arrow.
2. Insert the two AA rechargeable batteries. Be sure to match the positive (+) and negative (-) ends to those indicated on the battery cavity.
3. Slide battery compartment cover back into locked position.
4. Charge batteries following Battery Charging instructions.
5. Insert the headphone/earbud connector into the jack on the left side of the amplifier.
6. Be sure the volume control is in the "low" position before putting on the headphones/earbuds and adjusting for comfort.
7. Slowly, turn the On/Off/Volume control wheel. You will feel a "click" and the green LED light will illuminate when the SuperEar is turned on.
8. Continue to turn the volume control knob to the desired level of amplification.
9. Adjust the Frequency Selector to Low, Medium or High to filter the undesirable frequency ranges.
10. This amplifier has a maximum gain of over 50 dB. At approximately 1/4 open on the volume control, there will be about a 10dB sound gain. The maximum gain of well over 50 dB is delivered with the volume wheel all the way open. You may not need full volume gain. Choose an amplification level that is most comfortable for you.
11. If headphones/earbuds not in use are placed near the amplifier while it is turned on, feedback squeal will result.

Instructions for Use

1. Remove the battery compartment cover by lightly depressing the cover with your thumb and slide it in the direction of the arrow.
2. Insert the two AA rechargeable batteries. Be sure to match the positive (+) and negative (-) ends to those indicated on the battery cavity.
3. Slide battery compartment cover back into locked position.
4. Charge batteries following Battery Charging instructions.
5. Insert the headphone/earbud connector into the jack on the left side of the amplifier.
6. Be sure the volume control is in the "low" position before putting on the headphones/earbuds and adjusting for comfort.
7. Slowly, turn the On/Off/Volume control wheel. You will feel a "click" and the green LED light will illuminate when the SuperEar is turned on.
8. Continue to turn the volume control knob to the desired level of amplification.
9. Adjust the Frequency Selector to Low, Medium or High to filter the undesirable frequency ranges.
10. This amplifier has a maximum gain of over 50 dB. At approximately 1/4 open on the volume control, there will be about a 10dB sound gain. The maximum gain of well over 50 dB is delivered with the volume wheel all the way open. You may not need full volume gain. Choose an amplification level that is most comfortable for you.
11. If headphones/earbuds not in use are placed near the amplifier while it is turned on, feedback squeal will result.

Belt Clip Installation/Removal

- 1) Place belt clip over the Battery Compartment Door, with the top hinge centered between the two anti-scurf pads and the clips over the retaining slots.
- 2) Slip the small joint into the slot located above the Battery Compartment Door.
- 3) Apply pressure on both sides of the clip until the clip retainers snap into the slots on the sides of the unit.
- 4) To remove, flex out one side of the clip retainer, unsnap from the retaining slot and remove the clip.

Battery Charging

To charge the batteries, place the rechargeable batteries in the battery compartment. Plug the recharging Mini USB end of the Recharging cable into the recharging port. The other end of the cable with the USB connector may be either plugged into a USB port on your computer or plugged into the Recharging Wall Plug, which is in turn plugged into a 110VAC outlet. Once charging begins, the red Recharging Indicator LED will flash until the LED turns to solid red. This indicates charging is complete.

- Allow Ni-MH Batteries (provided) to charge 12-16 hours for up to 40 hours useage.
- Excessive charging may reduce the life expectancy and capacity
- Normal operation may continue while recharging, if the "Power" mode is ON.

Belt Clip Installation/Removal

- 1) Place belt clip over the Battery Compartment Door, with the top hinge centered between the two anti-scurf pads and the clips over the retaining slots.
- 2) Slip the small joint into the slot located above the Battery Compartment Door.
- 3) Apply pressure on both sides of the clip until the clip retainers snap into the slots on the sides of the unit.
- 4) To remove, flex out one side of the clip retainer, unsnap from the retaining slot and remove the clip.

Battery Charging

To charge the batteries, place the rechargeable batteries in the battery compartment. Plug the recharging Mini USB end of the Recharging cable into the recharging port. The other end of the cable with the USB connector may be either plugged into a USB port on your computer or plugged into the Recharging Wall Plug, which is in turn plugged into a 110VAC outlet. The charging cube is 110/240VAC capable and can be used with a simple plug tang adapter for any non-US market. Once charging begins, the red Recharging Indicator LED will flash until the LED turns to solid red. This indicates charging is complete.

- Allow Ni-MH Batteries (provided) to charge 12-16 hours for up to 40 hours useage.
- Excessive charging may reduce the life expectancy and capacity
- Normal operation may continue while recharging, if the "Power" mode is ON.

**This Device passed the FCC Verification Testing.
The unique FCC Identifier for this device is SE9000.**

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Specifications

- Designed to provide up to 50+ decibels of ambient sound amplification. Maximum output of 107dB.
- Frequency response: 100 to 14000 hertz.
- Solid state microelectronic circuitry requires only two Rechargeable AA (1.2V) batteries.
- Premium stereo earphones and earbuds deliver clear, crisp sound particularly well in the human voice range.
- Operating temperature range 0°C min, 55°C max.
- Relative humidity operating range: 20°C/95%RH

3 YEAR WARRANTY

Sonic Technology Products, Inc. warrants to the original purchaser that your SuperEar® will be free of defects in materials or workmanship for a period of 3 years from the date of purchase. Warranty is void if defects are caused by negligence, abuse, misuse or damage caused by water or heat. Returns made under this Warranty must be accompanied by \$3.00 to cover the handling and return postage cost. This Warranty gives you specific legal rights and you may have other rights which may vary from state to state. Should it be necessary to send your SuperEar® in for service under the terms of this Warranty, return it in a well-packaged carton and ship postage prepaid.

SonicTechnology
P.O. Box 539
Grass Valley, CA 95945
www.sonictechnology.com
11/2015

**This Device passed the FCC Verification Testing.
The unique FCC Identifier for this device is SE9000.**

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Specifications

- Designed to provide up to 50+ decibels of ambient sound amplification. Maximum output of 107dB.
- Frequency response: 100 to 14000 hertz.
- Solid state microelectronic circuitry requires only two Rechargeable AA (1.2V) batteries.
- Premium stereo earphones and earbuds deliver clear, crisp sound particularly well in the human voice range.
- Operating temperature range 0°C min, 55°C max.
- Relative humidity operating range: 20°C/95%RH

3 YEAR WARRANTY

Sonic Technology Products, Inc. warrants to the original purchaser that your SuperEar® will be free of defects in materials or workmanship for a period of 3 years from the date of purchase. Warranty is void if defects are caused by negligence, abuse, misuse or damage caused by water or heat. Returns made under this Warranty must be accompanied by \$3.00 to cover the handling and return postage cost. This Warranty gives you specific legal rights and you may have other rights which may vary from state to state. Should it be necessary to send your SuperEar® in for service under the terms of this Warranty, return it in a well-packaged carton and ship postage prepaid.

SonicTechnology
P.O. Box 539
Grass Valley, CA 95945
www.sonictechnology.com
11/2015

SuperEar®
Model SE9000
Personal Sound Amplifier



Owner's Manual

SuperEar®
Model SE9000
Personal Sound Amplifier



Owner's Manual